## MEMORANDUM

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### NT-PRO-BNP Collections Instructions – Central Laboratory, Glen Site

Please note that for the analysis of NT-proBNP at the Central Laboratory – Glen Site, which will go-live on Tuesday, February 22, 2022, the preferred tube type is a mint green tube (lithium heparin with gel).

As a reminder, IMMEDIATELY after the tube is pulled off from the holder, it must be gently inverted to ensure proper mixing. Proper mixing prevents plasma and/or serum to contain fibrin which could be aspirated by the analyzer(s) and cause erroneous results, especially for patients receiving anticoagulant or thrombolytic therapy (incomplete clot formation).

Transferring blood from any syringe into a vacutainer can also cause fibrin and is strongly discouraged.

Additionally, to limit hemolysis, butterfly needles should only be used when it is difficult to procure blood from a patient with a larger needle.

As of February 22th, 2022, the new orderable for Oasis users will change to NT PRO BRAIN NATRIURETIC PEPTIDE.