NOTE DE SERVICE MEMORANDUM

Date 2025-05-13

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To: **All Optilab-MUHC Laboratory Users**

De/From: **Dr Raymond Tellier MD MSc FRCPC CSPQ D(ABMM) Dr Gerasimos J. Zaharatos, MD, FRCPC(C)**

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Title Site Chief Microbiology Laboratory MUHC Glen Site Chief of Medical Microbiology, OPTILAB – Montréal MUHC

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Subject: **New kits and platform for cytomegalovirus (CMV) quantitative viral load testing**

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Dear colleagues,

Please note that **as of May 6, 2025** the MUHC (Glen Site) Microbiology laboratory will be transferring the CMV viral load assay from the Abbott Real-time CMV assay to a new commercial kit and platform: *Altostar* CMV 1.5 CAN (Altona Diagnostics). The new procedure has been extensively verified and correlates very well with the assay currently in use (but please note specific characteristics of the new assay below). This change is necessary because of the phasing out of the assay currently in use.

The *Altostar* CMV 1.5 CAN diagnostic assay is Health Canada - approved, and provide quantitative results expressed using the International Unit (IU) standard for CMV.

The *Altostar* CMV 1.5 CAN assay for expresses results in IU/mL; our validation has established an approximate conversion factor between IU/mL and the copies/mL results provided by Abbott assay as follows: **1.0 copies/mL (Abbott) = 6.60 IU/mL (Altostar)**

To give some examples:

- A result obtained with the Abbott assay of 1200 copies/mL (3.08 log copies/mL) would be approximately 7920 IU/mL (3.90 log IU/mL) with the *Altostar* assay
- A result of 67000 IU/mL (4.83 log IU/mL) with the *Altostar* would have been approximately 10152 copies/mL (4.01 log copies/mL) with the Abbott Assay

Please refer to the sample collection manual for information regarding sample requirements and transportation; there are no changes to these requirements with the new assay, which is done on plasma obtained from EDTA anticoagulated blood samples.

For any concerns or question regarding these new tests, please **contact Dr. Raymond Tellier** at the MUHC (Glen site) Microbiology Laboratory.

Thank you for your understanding