NOTE DE SERVICE

MEMORANDUM

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
Optilab Mtl - MUHC, Glen Site wishes to remind you that tubes drawn for coagulation testing must respect the anti-coagulant to blood ratio ($\pm 10\%$).

Recently, the laboratory has received an increase number of under-filled tubes. Upon further inspection, it was noted that the tubes used were expired which is the likely cause for loss of vacuum leading to under-filled tubes. **Please verify your stock and ensure that you do not use expired procurement tubes.**

On Sodium Citrate tubes, manufactured by BD Life Sciences, the minimum draw line (90% filled) is denoted by an etched line (minimum fill indicator) around the tube. To reach minimum level, the blood must be above or at the bottom edge of the minimum fill indicator.

If you are using a butterfly/winged collection set you must draw a discard tube to remove air from the collection line before drawing any tubes for coagulation studies.

If you are using a syringe to fill tubes be sure not to overfill the tubes. On a 2.7 mL tube the overfilled mark is denoted by the bottom edge of the cap, on a 1.8 mL tube the overfilled mark is 3 mm below the bottom edge of the cap (please refer to below diagram).

Any coagulation tubes that are either underfilled or overfilled will be cancelled as inadequate specimen.

Should you have any questions regarding this matter, please contact Client Services for Clinical Laboratories at 514-934-1934 ext 35687.

Thank you for your cooperation and understanding.