NOTE DE SERVICE

MEMORANDUM

Date 2024-11-11

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see below for English version

Dear colleagues,

Effective since January 1, 2025, to increase our sample management efficiency and reduce pre-analytical risks:

- 1) FFPE blocks will no longer be accepted at the MUHC Core Molecular Diagnostics Laboratory - Glen site. Only samples cut in advance, in accordance with our procedures, will be accepted.
- 2) H&E slides will not be returned as stipulated in our request.

To respect the conformity of the FFIP sections of molecular tests, it is imperative to work in such a way as not to contaminate the samples between them (cross-contamination).

The conformity of the cuts between each sample must follow the instructions:

The microtome blade must be cleaned with an RNase AWAY type product. Then we proceed to cut:

- BRCA: 10 scrolls with a thickness of 10 microns.
<https://cusm.ca/sites/default/files/docs/m-Labs/brca-somatic-instruction-test-description-fr.pdf>
- RNA-seq: 10 scrolls with a thickness of 10 microns <https://muhc.ca/sites/default/files/docs/m-Labs/targeted-rnaseq-instruction-test-description-en.pdf>
- Ampliseq Focus Panel: 6 scrolls with a thickness of 5 microns <https://muhc.ca/sites/default/files/docs/m-Labs/ampliseq-focus-panel-test-information-instructions-en.pdf>

When the sample has been placed in an Eppendorf tube, we clean the microtome with the RNase Away. We can then proceed for a second patient.

The information sheet for each test as well as the procedure for cutting the paraffin scrolls can be found in the link below and it is also attached to this memo.

<https://muhc.ca/sites/default/files/docs/m-Labs/ampliseq-focus-panel-test-information-instructions-en.pdf>

<https://muhc.ca/sites/default/files/docs/m-Labs/brca-somatic-instruction-test-description-en.pdf>

<https://muhc.ca/sites/default/files/docs/m-Labs/targeted-rnaseq-instruction-test-description-en.pdf>

For any questions please contact the Molecular Diagnostics laboratory of the MUHC Site Glen at 514-934-1934 x24900

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AmpliSeq Focus Panel – Preparation and Shipping Guidelines

The quality of the *AmpliSeq Focus Panel* result depends on the following **critical pre-analytical steps** in order to obtain high-quality nucleic acids. This includes fixation, microtome section preparation, specimen preparation and shipping.

1. Fixation:

Fixation of the specimen with only **10% neutral-buffered formalin**, for a period of 24–72 hours, allows preservation of nucleic acids. Other fixatives may lead to false negative results and are **not** accepted (e.g. Bouin, B5, AZF, Holland's or strong acid decalcification procedures). Under fixation or over fixation can lead to suboptimal results.

2. Microtome Section Preparation:

Ribonucleases (RNases), which is present ubiquitously (eg. normal human skin), can destroy nucleic acids. Careful preparation and clean material are necessary to avoid cross contamination and RNases. Disposable material and smooth tip tweezers cleaned with an ANTI-RNase product (eg. RNase Away™) are encouraged.

- 1) Wear gloves and change them frequently, especially between specimens
- 2) Change the microtome blade for each case
- 3) Use clean instruments to transfer the sections (i.e. sterile needle to transfer the scrolls into clean 1.5 ml Eppendorf tubes)
- 4) Use hot plate or a clean water bath (for unstained slides)
- 5) All material (tube, slides) must be adequately labelled (two patient identifiers)

3. Specimen Preparation:

Specimen type	Requirements
Scrolls FFPE*	<p>≥ 10% of tumor cell content (TCC) in the whole cut surface is required. TCC is the number of viable nucleated tumor cells in proportion to the total number of cells.</p> <p>1) 10 scrolls at 5-µm (minimum 5 scrolls) in a labelled 1.5 ml Eppendorf tube</p> <p>2) 1 H&E slide</p>
Unstained Slides FFPE	<p>Required for microdissection when TCC <10%</p> <p>1) 6 unstained labelled slides at 5-µm thickness, uncharged and unbaked</p> <p>2) 1 H&E slide, pointed for tumor area if tumor content < 10%</p>

*FFPE: Formalin-fixed paraffin-embedded. When the tumor cell content (TCC) ≥ 10%, scrolls are preferred.

AmpliSeq Focus Panel – Test Information

Indications for ordering

This test allows the detection of genomic abnormalities with a diagnosis, prognosis or predictive role in diverse solid tumours (e.g. lung, lower GI, melanoma).

Test Methodology

The *AmpliSeq Focus Panel* is a targeted next generation sequencing assay for analysis of hotspot regions in 52 genes with known relevance to solid tumours. The panel enables the simultaneous analysis of genes associated with several most common cancers (e.g. lung, colon). The DNA and RNA analysis can detect multiple alterations including: (i) single nucleotide variants (SNVs), (ii) small insertion / deletions (indels), (iii) copy number gains and (iv) gene fusions in RNA samples. See below for gene list. For more information on the hotspot regions tested, please contact the laboratory.

Gene List

DNA (hotspot variations and *amplifications)				
AKT1	EGFR*	FGFR4*	JAK3	MYCN*
ALK*	ERBB2*	GNA11	KIT*	NRAS
AR*	ERBB3	GNAQ	KRAS*	PDGFRA*
BRAF*	ERBB4	HRAS	MAP2K1	PIK3CA*
CCND1*	ESR1	IDH1	MAP2K2	RAF1
CDK4*	FGFR1*	IDH2	MET*	RET
CDK6*	FGFR2*	JAK1	MTOR	ROS1
CTNNB1	FGFR3*	JAK2	MYC*	SMO
DDR2				

RNA (fusions)				
ABL1	EGFR	ETV5	NTRK1	PPARG
ALK	ERBB2	FGFR1	NTRK2	RAF1
AKT3	ERG	FGFR2	NTRK3	RET
AXL	ETV1	FGFR3	PDGFRA	ROS1
BRAF	ETV4	MET		

Test Interpretation

Only variants with a recognized clinical significance, diagnostic, prognostic or predictive are reported (Tier I and II; PMID: 27993330). Additional results are available upon request.

Note: This test does not detect large deletions.

Type of specimens accepted

- Tumour cell content $\geq 10\%$ is required. This information is mandatory to assess the validity of the test.
- Specimen (cytology or histology), frozen or formalin-fixed paraffin-embedded (FFPE) are accepted.
- For FFPE: only SCROLLS or UNSTAINED SLIDES are accepted.
- Unstained slides: a minimum of 5 sections of 5 μm thickness non-baked and uncharged, with 1 H&E slide (tumour area should be outlined with permanent marker if macrodissection is needed).

The full procedure is detailed in the document Preparation and Shipping Guidelines.

Limitations

Results must be interpreted in the **context of clinical, radiological and histological findings**. If results obtained do not match other clinical or laboratory findings, or if you have novel relevant information, please contact the laboratory as soon as possible for updated interpretation.

Results may be compromised if the recommended procedures (tissue fixation and preparation) have not been followed.

A negative (wild type) result does not fully rule out the presence of an alteration but may be linked with limits of detection of this assay (i.e. insufficient % of tumour cell content or poor fixation).

Rare polymorphisms may be present that could lead to false-negative or false-positive results.

Turnaround time: 10 working days

Revised & approved by Dr. Gomez 23-10-26

BRCA Tumour Panel – Sample Preparation and Shipping Guidelines

The quality of the *BRCA Tumour Panel* result depends on the following **critical pre-analytical steps** in order to obtain high-quality nucleic acids. This includes fixation, microtome section preparation, specimen preparation and shipping.

1. Fixation:

Fixation of the specimen with only **10% neutral-buffered formalin**, for a period of 24–72 hours, allows preservation of nucleic acids. Other fixatives may lead to false negative results and are not accepted (e.g. Bouin, B5, AZF, Holland's or strong acid decalcification procedures). Under fixation or over fixation can lead to suboptimal results.

2. Microtome Section Preparation:

Ribonucleases (RNases), which is present ubiquitously (e.g. normal human skin), can destroy nucleic acids. Careful preparation and clean material are necessary to avoid cross contamination and RNases. Disposable material and smooth tip tweezers cleaned with an ANTI-RNase product (eg. RNase Away™) are encouraged.

- Wear gloves and change them frequently, especially between specimens
- Change the microtome blade for each case
- Use clean instruments to transfer the sections (i.e. sterile needle to transfer the scrolls into clean 1.5 ml Eppendorf tubes)
- Use hot plate or a clean water bath (for unstained slides)
- All material (tube, slides) must be adequately labelled (two patient identifiers)

3. Specimen Preparation:

<i>Specimen type</i>	<i>Requirements</i>
Scrolls FFPE*	TCC is the number of viable nucleated tumour cells in proportion to the total number of cells. <ul style="list-style-type: none"> • ≥ 10% of tumour cell content (TCC) in the whole cut surface is required. • 10 scrolls at 10-µm (minimum 5 scrolls) in a 1.5 ml Eppendorf tube with 2 patients' identifiers should be submitted • 1 H&E slide
Unstained Slides FFPE*	Required for macrodissection when overall TCC <10% <ol style="list-style-type: none"> 1) 6 unstained labelled slides at 10-µm thickness, uncharged and unbaked 2) 1 H&E slide, pointed for tumour area if TCC < 10%

(*) FFPE: Formalin-fixed paraffin-embedded.

Revised 23-10-26

BRCA Tumour Panel – Test Information

Indications for ordering

This test detects genomic abnormalities in the *BRCA1* and *BRCA2* genes that may give the patient eligibility for targeted therapy with PARP inhibitors. This test may be ordered for patients with the following tumour types:

- Breast (Triple-negative, HER2-negative, high-risk ER+, gBRCA-positive*)
- Ovarian (High-grade non-mucinous)
- Prostate (Castration-resistant)
- Pancreas (all)

*Patients with a germline *BRCA1/2* pathogenic variant

Test Methodology

The *BRCA Tumour Panel* is a targeted next generation sequencing capture panel for analysis of the entire *BRCA1* and *BRCA2* genes. This test only analyzes DNA and can detect multiple types of alterations including: (i) single nucleotide variants (SNVs), (ii) small insertion/deletions (indels), (iii) copy number variants (deletions and duplications).

Genes Tested

- *BRCA1*
- *BRCA2*

Test Interpretation

Only variants classified as Pathogenic or Likely pathogenic according to modified ACMG criteria are reported (PMID: 25741868, <https://cspec.genome.network/cspec/ui/svi/affiliation/50087>). Additional results are available upon request. This test cannot definitively determine whether a variant is present in the germline or restricted to the tumour.

Specimens accepted

- **Tumour cellularity $\geq 10\%$** is required. This information is mandatory to assess the validity of the test.
- **Scrolls or unstained slides** from Formalin-fixed paraffin-embedded (FFPE) specimens (cytology or histology) are accepted
- **FFPE blocks** are NOT accepted
- **Scrolls: Ten scrolls at 10um**
- **Unstained slides:** 5 unstained labelled slides of 10 μm thickness, unbaked and uncharged, with 1 H&E slide (tumour area should be outlined with permanent marker if macrodissection is needed).
- The full procedure is detailed in the document "*Sample Preparation and Shipping Guidelines*".

Limitations

Results must be interpreted in the context of clinical, radiological and histological findings. If results obtained do not match other clinical or laboratory findings, or if you have novel relevant information, please contact the laboratory as soon as possible for updated interpretation.

Results may be compromised if the recommended procedures (tissue fixation and preparation) have not been followed. A negative result does not fully rule out the presence of an alteration and may be due to the limitations of this assay (i.e. insufficient % of tumour cell content or poor fixation).

CNV detection cannot be guaranteed if the tumour cellularity is $<40\%$

Turnaround time: 4-6 weeks

Targeted RNASeq Panel – Test Information

Indications for ordering

This test detects translocations/fusions, exon skipping, internal tandem duplications, and sequence variants in select oncogenes (see **Gene list**) that may have diagnostic, prognostic, or therapeutic implications for patients with certain tumour types (see below). This test may be ordered for patients with the following tumour types:

- Sarcomas
- Cholangiocarcinomas
- Brain tumours
- Rare solid tumours
- Certain hematological malignancies

Test Methodology

The *Targeted RNASeq panel* is a targeted next generation sequencing capture panel for analysis of a targeted list of genes, relevant to certain tumour types. This test only analyzes RNA and can detect multiple types of alterations, including (i) translocations/fusions, (ii) internal tandem duplications, (iii) single nucleotide variants (SNVs) and small insertion/deletions (indels), (iv) alterations that cause exon skipping.

Test Interpretation

Only variants with a recognized or predictive diagnostic, prognostic or therapeutic clinical significance are reported (Tier I and II; PMID: 27993330). Additional results are available upon request. This test cannot definitively determine whether a variant is present in the germline or restricted to the tumour.

Specimens accepted

- **Tumour cellularity $\geq 20\%$** is required. This information is mandatory to assess the validity of the test.

Preferentially:

- **Frozen tissue** minimally 30 mg of tissue on dry ice
- **Fresh tissue** on RNA later minimally 30 mg of tissue

If no fresh/frozen tissue is available:

- **FFPE blocks** are NOT accepted
- **Scrolls or unstained slides** from Formalin-fixed paraffin-embedded (FFPE) specimens (cytology or histology) are accepted

Scrolls: Ten scrolls at 10um

Unstained slides: 10 unstained labelled slides of 10 μm thickness, unbaked and uncharged, with 1 H&E slide (tumour area should be outlined with permanent marker if macrodissection is needed).

- The full procedure is detailed in the document "*Sample Preparation and Shipping Guidelines*".

Limitations

Results must be interpreted in the context of clinical, radiological and histological findings. If results obtained do not match other clinical or laboratory findings, or if you have novel relevant information, please contact the laboratory as soon as possible for updated interpretation.

Results may be compromised if the recommended procedures (tissue fixation and preparation) have not been followed.

A negative result does not fully rule out the presence of an alteration and may be due to the limitations of this assay (i.e. insufficient % of tumour cell content or poor fixation).

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