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<b>Sector:</b> Anatomic Pathology	
<b>Topic:</b> Cytology	
<b>Location :</b> <a href="#">Site Web d'OMNI au CUSM</a>	<b>Effective date:</b> 2024-11-20 <b>Review Period:</b> 3 years
<b>Key words:</b> specimen collection, cytology	<b>Distribution:</b> <input checked="" type="checkbox"/> All pathology laboratory personnel and clients
<b>Version number:</b> 1.0	<b>Annex:</b>
<b>Written by / Signature:</b> Dr. Manon Auger Julie Brodeur Bianca Blais-Bourgoin	<b>Date (yyyy/mm/dd):</b> 2024-07-24 2024-07-24 2024-07-24
<input type="checkbox"/> Electronically signed document	
<b>Reviewed by / Signature:</b> Jennifer Pamplin	<b>Date (yyyy/mm/dd):</b> 2024/10/16
<input type="checkbox"/> Electronically signed document	
<b>Authorized by / Signature(s):</b> Melissa Trickey	<b>Date (yyyy/mm/dd):</b> 2024/11/06
<input type="checkbox"/> Electronically signed document	
<b>Replaces the version:</b> New document	
<b>Changes to the last authorized version:</b>	

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**Note:** This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.

Date effective: 2024-11-20	Code ID: SOP-01-APHI-PC-001075 EN	Version: 1.0
----------------------------	-----------------------------------	--------------

1 Application .....	2
2 Objective .....	2
3 Definitions .....	2
4 References.....	2
5 Responsibilities .....	2
6 Equipment and materials .....	2
7 Procedure .....	2
8 Related documents .....	12

## 1 Application

This procedure provides instructions for the collection of delivery of cytological samples to the laboratory.

## 2 Objective

This procedure aims to describe the steps necessary to provide a sample to the cytopathology laboratory.

## 3 Definitions

## 4 References

## 5 Responsibilities

- Management is responsible for providing appropriate information to the client and keeping all policies and procedure up to date.
- All clients are responsible for knowing and following the contents of this protocol.

## 6 Equipment and materials

## 7 Procedure

### **LABORATORY HOURS, SPECIMEN COLLECTION, PHONE NUMBERS AND LOCATION**

#### **CYTOPATHOLOGY LABORATORY HOURS:**

The Cytopathology Laboratory is open to receive specimens from 8:00 am until 16:00 pm, Monday through Friday. It is closed on weekends and statutory holidays recognized by the MUHC.

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

**SPECIMEN COLLECTION:****RVH SITE:**

All cytology specimens should be sent to the lab using the MUHC/Royal Victoria Hospital Transport service. The lab receives specimens, via the porter service, four (4) times per day. If the specimen is deemed to be STAT, do not put it out for routine transport (consider sending it with an orderly or bringing it yourself).

**MGH SITE:****Cytology Specimens can be dropped off**

All specimens for Cytology are to be deposited in the refrigerator located in the Pathology Department room C3-152.

Non-gyn specimens (consisting of no more than 100ml of any liquid) are to be placed in orange-topped specimen bottles (with 15 to 30 ml of Cytolyt please), the caps well-sealed, and the jars placed in plastic bags (to contain any leakage). The bottles and requisitions are to be well identified as to patient, specimen origin, treating staff physician, location (ward, clinic, office).

PLEASE NOTE: Any large quantities of liquid received (over and above the 100 ml requested) will be returned to the ward/clinic for downsizing.

If it is necessary to collect a fluid after working hours, place it in a refrigerator overnight.

**SEND THE SPECIMEN TO THE CYTOLOGY LABORATORY EARLY IN THE MORNING OF THE NEXT WORKING DAY.**

Turn Around Time (TAT): NON-GYN: Within 2-3 working days

GYN – Colposcopy and high risk patient: Within 3-5 working days  
GYN – Routine: Variable

**PHONE NUMBERS:**

**External: 514-934-1934 ext: 38782**

**514-934-1934 ext: 35687 option 2 (for results request)**

**Internal: 22803 (reports)**

**38784 (prep lab, supplies, specimen info)**

**38782 (Cytology Co-Coordinator – all info)**

**LOCATION OF THE CYTOLOGY LABORATORY:**

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

## 7.1 GYNECOLOGICAL CYTOLOGY

### GENERAL PRINCIPLES IN OBTAINING SPECIMENS:

1. Instruct the patient not to use a vaginal douche for 24 hours prior to appointment (may wash away diagnostic cells).
2. Let at least four (4) weeks lapse post surgery (cryotherapy, cone biopsy, LASER, D & C) and six (6) weeks post radiotherapy, to allow for sufficient healing. Smears taken too early may show marked reparative changes that are difficult to differentiate from malignant changes.
3. In general, allow 4-6 weeks to elapse after an abnormal smear, before doing a repeat.
4. Moisten the speculum with warm tap water only. Lubricants distort and obscure the cells in a stained smear.
5. Take the smear before performing a manual examination. The exam could dislodge the cells, contaminate the smear with lubricant and/or glove powder.

**NOTE:** Dipping the spatula into saline first will help prevent air-drying artefact in those women with dry vaginas or when taking vulvar smears

### MATERIALS REQUIRED:

1. Clean glass slide with one frosted end.
2. Cytology spray fixative.
3. Speculum
4. Wooden cervical spatula (Duo-tip)
  - Cervical cytobrush (optional)
5. Cardboard slide holder or slide box
6. Cytology requisition form (DM-5048 or DM-5049)

### METHODS OF OBTAINING CONVENTIONAL SMEARS

#### Cervical Smear:

1. Label slide on frosted end of glass in pencil, **with two patient identifiers**. Complete name (first and last name) **and** date of birth or medical record number.

#### ***DOUBLE IDENTIFICATION IS REQUIRED***

2. Place fixative in an easily accessible location.
3. Insert speculum, using water as a lubricant if necessary.
4. Place pointed end of Duo-top spatula as high into the cervical canal as possible, to ensure contact with the transformation zone.
5. Press firmly and rotate the spatula 360<sup>0</sup> to scrape the entire squamo-columnar junction.
6. Transfer the material to a labeled glass slide, spreading evenly.
7. **FIX IMMEDIATELY** with spray fixative ( 2 pumps of the fixative are generally sufficient).
8. Send in box (if many slides are being sent) or cardboard mailer (if only the occasional slide is being sent) along with a properly and legibly completed Cytology requisition form

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

Endocervical Brush Smear:

1. Label slide on frosted end of glass in pencil, **with two patient identifiers** Complete name (first and last name) **and** date of birth or medical record number.

**DOUBLE IDENTIFICATION IS REQUIRED**

2. Insert brush gently into the endocervix and slowly rotate 90°.
3. Prepare a smear by rolling and twisting the brush on the slide
4. Immediately apply spray fixative.
5. Send to the cytopathology lab along with a completed cytology requisition.

Vaginal Smear:

1. Label slide on frosted end of glass in pencil, **with two patient identifiers** Complete name (first and last name) **and** date of birth or medical record number.

**DOUBLE IDENTIFICATION IS REQUIRED**

2. Insert speculum, using water as lubricant if necessary.
3. Gently scrape the upper third of the vaginal wall with a wooden cervical spatula or tongue depressor.
4. Spread material evenly on the glass slide.
5. Spray fix immediately.
6. Send to the cytopathology lab along with a completed cytology requisition.

Vulvar Smear:

1. Label slide on frosted end of glass in pencil, **with two patient identifiers** Complete name (first and last name) **and** date of birth or medical record number.

**DOUBLE IDENTIFICATION IS REQUIRED**

2. Dip the flat end of wooden spatula or tongue depressor in saline or water (to minimize air- drying artefact).
3. Scrape vulva.
4. Spread material evenly on the glass slide.
5. Spray fix immediately.
6. Send to the cytopathology lab along with a completed cytology requisition.

ThinPrep specimens (Cervix)

1. With a combined brush or spatula
  - a. Insert the brush in the endocervix canal until the short part of the brush is in contact with the endocervical zone.
  - b. Push gently and turn 5 times clockwise.
  - c. Rinse the brush 10 times in the PreservCyt solution.
  - d. Turn vigorously the brush to remove most of the cellular material. e- Put the cap
  - e. Put the cap
  - f. Label the container with the patient full name and identification number or RAMQ.

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**Note:** This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

***DOUBLE IDENTIFICATION IS REQUIRED.***

2. With endocervix brush
  - a. Insert the brush in the cervix.
  - b. Turn slowly  $\frac{1}{4}$  or  $\frac{1}{2}$  turn.
  - c. Rinse the brush 10 times in the PreservCyt solution.
  - d. Shake vigorously the brush to remove most of the cellular material.
  - e. Put the cap
  - f. Label the container with the patient full name and identification number or RAMQ.

***DOUBLE IDENTIFICATION IS REQUIRED.*****Anal Smear (Preference: use Thin Prep--- but conventional smears are also acceptable)**

1. Label slide on frosted end of glass in pencil, **with two patient identifiers** Complete name (first and last name) **and** date of birth or medical record number.

***DOUBLE IDENTIFICATION IS REQUIRED***

2. Moisten the swab (not cotton) with water. Do not use lubricant.
3. The swab should be inserted approximately 5 to 6 cm into the anal canal. Once inserted deep enough into the anus (necessary in order to collect both rectal columnar and anal squamous cells) the swab should be pulled out, applying some pressure to the wall of the anus, rotating the swab in a spiral motion along the way.
4. Spread on a glass slide.
5. Spray fix immediately.
6. Send to the cytopathology lab along with a completed cytology requisition.

**7.2 BODY CAVITY FLUID CYTOLOGY****GENERAL PRINCIPLES:**

1. 100 ml is the ideal volume of fluid to be sent to the Cytology lab ( as little as 2-3 ml can be processed if that is all that is available). **NEVER SEND MORE THAN 100 ml.** Do not send large vacuum 'tap' bottles. If necessary to remove liquid from these large containers, shake the jars to resuspend the cells, then aspirate out a maximum of 100 ml into a wide mouthed jar and send this to the lab.
2. Long -standing effusions frequently contain large numbers of inflammatory cells, often making it difficult to visualise cells and therefore render a diagnosis. A repeat tap, 24 hours after the first tap, is usually sufficient to greatly improve the "view".

**MATERIALS REQUIRED:**

1. Paracentesis tray.
2. Clean, dry, non-sterile wide mouthed jar with tight fitting lids
3. Cytology requisition form (DM-5047)

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

**METHOD OF OBTAINING SPECIMEN:**

1. Perform para/thoracentesis. If patient is bedridden, it is advisable to gently rotate him prior to tapping to re-suspend cells in the body cavity.
2. Collect up to 100 ml of body cavity fluid and place into a clean, non-sterile wide-mouthed jar. **DO NOT ADD FIXATIVE TO THESE FLUIDS.**
3. Immediately send specimen accompanied by a completed cytology requisition, to the cytology lab. Place the specimen in the refrigerator if there is going to be any delay in transport.

**7.3 RESPIRATORY CYTOLOGY****TYPES OF SPECIMENS:**

1. Sputum series
2. Post-bronchoscopy sputum
3. Aerosol-induced sputum
4. Bronchoscopy specimens
5. Tracheal Aspirates
6. Fine Needle Aspiration of the Lung

**GENERAL PRINCIPLES IN OBTAINING SPECIMENS:**

1. A series of 3-5 early morning deep cough sputa is recommended. The chance of detecting malignancy increases with the number of specimens examined.
2. Aerosol-induced sputa are indicated for patients who are incapable of unassisted sputum production.
3. Post-bronchoscopy sputum specimen is the single most important sputum. Obtaining this specimen should be a routine part of the bronchoscopy procedure. By definition, the post- bronchoscopy sputum is the first sputum produced immediately after a bronchoscopy has been performed.
4. Fine needle aspiration cytology of the lung should be considered in all patients whose pulmonary lesion remains obscure after routine investigations.

**METHODS OF OBTAINING SPECIMENS:****SPUTUM SERIES**

1. A series consists of 3-5 satisfactory, fresh, unfixed (if possible), early morning deep cough specimens.
2. **NEVER send** 24 hour sputum collections.

**MATERIALS REQUIRED:**

1. Wide mouth plastic jar (clean, not sterile) containing no fixative. However if the specimen is being taken on an out-patient basis, add 30 ml of Cytolyt, is to be added to the jar as a fixative.
2. Instruction sheet to be given to the patient.
3. Cytology requisition form (DM-5047)

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

**METHOD OF SPUTUM COLLECTION:**

1. The night prior to collecting a sputum series, give the patient a sputum cup and collection sheet.
2. Send specimen along with completed cytology requisition form to the cytopathology lab.
3. Repeat this procedure until 3-5 satisfactory specimens have been obtained.

**POST-BRONCHOSCOPY SPUTUM**

1. Give the patient a wide-mouthed jar before the bronchoscope is withdrawn.
2. Instruct the patient to cough deeply and to expectorate all sputum into the cup for the next hour.
3. Collect the cup after one hour, label the requisition form as “post-bronchoscopy sputum” after filling in the remainder of the information and send it and the specimen immediately to the cytology lab. For specimens that may encounter delays of more than 24 hours, 30 ml of Cytolyt or when unavailable and equal volume of 50% ethanol may be added to the jar.

**AEROSOL-INDUCED SPUTUM**

Arrangements for obtaining aerosol-induced sputa may be made by phoning Respiratory Therapy.

**BRONCHOSCOPY SPECIMENS**

1. Cytology specimens may be obtained during bronchoscopy by aspiration of secretions after BRONCHIAL WASHING with physiologic saline, and/or BRONCHIAL BRUSHINGS.
2. Specimens are placed in jars containing 30 ml Cytolyt, when possible, otherwise the specimens are sent to the lab, unfixed.
3. All specimens obtained should be labeled as to the patient and the site of procurement
4. Each container should be accompanied by its own requisition properly and legibly filled out.

**BRONCHOSCOPY SPECIMENS FOR PNEUMOCYSTIS JIROVECI (PJP)**

The following is the policy of the Cytopathology Laboratory for processing specimens for Pneumocystis jiroveci (PJP).

**For same day results:** The specimen must be **received** in the cytology lab **NO LATER THAN 2 PM (14:00hrs)**. If the specimen is received later than 2 p.m. results will only be available the next working day.

**7.4 URINARY TRACT CYTOLOGY****TYPES OF SPECIMENS:**

1. Voided urine specimens
2. Catheterized bladder specimens
3. Ureteral catheterization specimens
4. Specimens obtained at cystoscopy

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

### GENERAL PRINCIPLES OF URINARY TRACT CYTOLOGY:

1. Urinary tract cytology is primarily indicated where malignancies of urothelium are suspected. Malignancies of the kidney or prostate are rarely detected in routine urine specimens.
2. Cell morphology is better preserved in dilute rather than concentrated urine; so hydration of patients prior to obtaining specimen is encouraged, whenever possible. Cell degeneration occurs rapidly in urines. All specimens must be sent to the cytopathology laboratory AS SOON AS POSSIBLE. If a delay is anticipated, an equal volume of 50% ethyl alcohol may be added to the specimen. An equal volume of Cytolyt may also be added. (This is available at various locations in the MUHC hospitals – please call the lab (38782) for locations at the various sites).
3. Bladder and ureteral catheterization, or especially nephrostomy tubes, may induce cellular alterations which affect interpretation. Please fill in the requisition to indicate whether the specimen is voided or is obtained by catheterization
4. **NEVER SEND 24 HOUR URINE COLLECTIONS OR FIRST MORNING VOIDED URINES.**

### MATERIALS REQUIRED:

1. Clean dry container with tight fitting lid.
2. Appropriate urological instruments where indicated.
3. 50% Ethyl alcohol or Cytolyt, if needed
4. Cytology requisition form (DM-5047)

### TECHNIQUES OF SPECIMEN COLLECTION:

#### **VOIDED URINE**

1. Send the specimen fresh, if a minimal delay is anticipated. If a long delay (over 12 hours) is anticipated, fix voided urine in an equal volume of 50% ethyl alcohol (or Cytolyt). Please note on the the specimen which fixative is used..
2. Send this specimen immediately to the cytopathology lab along with a completed cytology requisition form.

#### **CATHETERIZED BLADDER SPECIMEN**

1. Send the specimen fresh, if a minimal delay is anticipated. If a long delay (over 12 hours) is anticipated, fix voided urine in an equal volume of 50% ethyl alcohol (or Cytolyt). Please note on the specimen which fixative is used.
2. Send this specimen immediately to the cytopathology lab along with a completed cytology requisition form.

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

## CATHETERIZED URETERAL SPECIMEN

1. Catheterize ureter to a point just below the suspected lesion.
2. Collect urine for ½ hour
3. Indicate on the collection bottle and the requisition form the side from which the specimen was obtained.
4. Send the specimen fresh, if a minimal delay is anticipated. If a long delay (over 12 hours) is anticipated, fix voided urine in an equal volume of 50% ethyl alcohol (or Cytolyt). Please note on the specimen which fixative is used.
5. Immediately send the specimen to the cytopathology laboratory along with a completed cytology requisition.

## CYTSOSCOPY SPECIMENS

1. Wash the bladder with physiological solution.
2. Send the washing along with a completed cytology requisition to the cytopathology lab.

## 7.5 CEREBROSPINAL FLUID CYTOLOGY (CSF)

### TYPES OF SPECIMENS:

1. Lumbar puncture fluid
2. Ventricular fluid

### GENERAL PRINCIPLES:

Cells degenerate very rapidly in CSF. Try to obtain during laboratory hours and arrange to have it sent to the cytology lab IMMEDIATELY after it has been obtained. If it is impossible to obtain the CSF during lab hours, refrigerate the specimen overnight and send to the lab as soon as it opens the next day.

### MATERIAL REQUIRED:

1. Spinal tap tray
2. Sterile screw cap jar or test tube.
3. Cytology requisition form (DM-5047)

### TECHNIQUE OF OBTAINING SPECIMEN:

1. Perform spinal tap
2. Collect as much fluid as possible into a sterile screw cap jar or test tube (even a few drops can be processed, but more is desirable).
3. Immediately send specimen and completed cytology requisition form to the laboratory.

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**Note:** This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

## 7.6 PROCEDURE FOR FINE NEEDLE ASPIRATION (FNA)

### 1. SUPERFICIAL/PALPABLE LESIONS

#### 1.1 PREPARATION

a. ANCILLARY STUDIES

Think about the possibility of need for ancillary studies, esp. microbiology cultures. The appropriate materials necessary for these ancillary studies should be immediately available for specimen handling after aspiration.

b. THE PATIENT

An appropriate history should be obtained and a directed physical examination should be performed as appropriate. The FNA technique and its potential complications should be clearly explained to the patient. Obtain the patient's verbal permission for the performance of the procedure.

c. THE MASS

Position the patient to ensure best access to the lesion, as well as the maximal safety and comfort of the patient. Palpate the mass carefully. It is often useful to palpate with the palmar surface of the four fingers of the dominant hand, using a circular motion. If the lesion cannot be localized, defer the study as a "blind" FNA has a low diagnostic yield. If the mass cannot be palpated but is radiologically demonstrable, then it is recommended to defer the study and proceed using radiologic or ultrasound guidance.

d. ANTISEPTICS

For most superficial FNAs, little preparation beyond handwashing, gloves and wiping the skin over the aspiration site with an alcohol swab is necessary. For the usual FNA, the preparation is comparable to that performed before inserting an intravenous line.

e. ANESTHESIA

Most superficial FNAs are performed without local anesthesia. In selected situations, or if multiple passes for a lesion are planned, superficial injection of a small amount of local anesthetic may be considered.

#### 1.2 EQUIPMENT

a. SYRINGES

The syringe that is chosen for this procedure should be chosen according to individual preference. There is no significant increase in the aspirating power of the larger syringe.

Before use, the syringe plunger should be moved to break any seal that might have occurred during manufacturing or storage.

b. NEEDLES

Disposable 22-gauge to 25-gauge needles with a clear hub are used. The clear hub makes it easier to see the aspirated material entering the hub. The length of the

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**Note:** *This is a document with a controlled distribution. It must be checked prior to use against the official version in the MUHC's OMNI website.*

Date effective: 2024-11-20

Code ID: SOP-01-APHI-PC-001075 EN

Version: 1.0

needles varies from 1.0 to 1.5 cm depending on the depth of the lesion.

c. **SYRINGE HOLDERS**

A number of different types of syringe holders are available. Their purpose is to enable the person performing the aspiration to comfortably direct the needle and draw back on the plunger of the syringe with only one hand, while stabilizing the mass to be aspirated with the other hand. Aspiration without a syringe holder is technically more difficult and can result in a less satisfactory specimen.

### **1.3 ASPIRATION TECHNIQUE**

a. **STABILIZE THE MASS**

Before puncture and during aspiration, the aspirator should stabilize the mass against deeper tissues with the index and middle fingers, or thumb and forefinger, depending on individual preference. To help stabilization, the skin over the lesion may be stretched.

b. **PIERCE THE SKIN**

Using a smooth motion, pierce the skin and advance the needle carefully into the mass, **WITHOUT APPLYING SUCTION**. A change in the texture of the tissue can be felt as a lesion is entered.

c. **APPLY SUCTION**

Apply and maintain negative pressure by drawing back on the plunger of the syringe. **DO NOT PUMP BACK AND FORTH ON THE SYRINGE PLUNGER DURING ASPIRATION.**

d. **SAMPLE LESION**

Move the needle back and forth in rapid 3- to 5-mm strokes. Discontinue aspiration after 10 to 20 strokes, or fewer if material is visible in the hub.

e. **NEEDLE WITHDRAWAL**

Release suction before removing the needle from the patient. Withdrawing the needle while maintaining suction will pull material back into the syringe, making the material virtually impossible to recover.

f. **COLLECTION OF MATERIAL:**

- i. Simply rinse thoroughly the content of the syringe and needle into a container containing Cytolyt (preferred), 50% ethyl alcohol or normal saline if transport to the lab is not delayed.
- ii. Send slides and container to the Cytology Laboratory with appropriate patient identification and an adequately filled requisition.

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