NOTE DE SERVICE

MEMORANDUM

date: 2026-06-15

à : Tous les clients des laboratoires de Microbiologie de la grappe Optilab-CUSM
to: *All clients of the Microbiology Laboratories of the Optilab-MUHC*

de: **Dr. Alan Spatz**
from: Directeur Médical/Medical Director
Optilab-CUSM/ Optilab-MUHC

Enzo Caprio
Directeur Clinico-Administratif/ Clinical-Administrative Director
Optilab-CUSM / Optilab-MUHC

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Dear Colleagues

As part of our preparedness for high-risk infectious diseases, including **Ebola virus disease (EVD)**, it is essential to reinforce strict requirements governing the handling of suspected biological specimens for pathogen from risk group 4 ([click here for more informations](#)).

In accordance with recommendations from the **INSPQ** ([cliquer ici](#)), all specimens suspected of Ebola must be managed under stringent biosafety protocols:

1. Safe Specimen Collection (at Bedside)

- Limit collections to **essential testing only**

Biochemistry	Hematology	Microbiology	Ebola (sendout)
Na, K, Cl, iCa, Anion Gap, Glu, BUN, Crea, Hct, TCO ₂ , Hgb	Manual cell count	Malaria (rapid test rapide and parasitemia) Hemoculture	LSPQ LNM

- **Blood collection tubes (in ordre):**
 - **1 tube SST** (yellow cap)
 - **1 lithium heparin or dosium tube** (dark green cap)
 - **4 EDTA tubes** (lavender cap)
 - **2 blood culture bottles.**

2. Transfusion medicine

- For a patient suspected of having EVD who requires a transfusion, the protocol for transfusion will be limited to O transfusion. No additional testing will be performed to determine the blood type until negative results for EVD are received.

3. Specimen Packaging

- Strict adherence to **triple packaging requirements** ([click here](#)):
 1. Leak-proof primary container (securely sealed)
 2. Leak-proof secondary container with absorbent material
 3. Rigid outer container properly labeled ([click here](#))
- Disinfection of external surfaces prior to transport
- Clear labeling:
 - **"Biohazard – Suspected Ebola"**
 - Proper patient identification

4. Safe Transport

Internal Transport (after disinfection of containers):

- **Prior notification** at the Laboratory before specimen arrival is mandatory
- Hand to hand delivery
- **No pneumatic tube systems usage**

External Transport:

- **Prior notification** at the Laboratory before specimen arrival is mandatory
- Compliance with Transportation of Dangerous Goods (TDG) regulations ([click here](#))
- Coordination with authorized carriers.
- Full documentation and ensured traceability.
- Transport of dangerous good (TDG) training is **MANDATORY** ([click here](#)).

Handling suspected Ebola specimens requires **absolute rigor at every step**.
Strict adherence to protocols is essential.