

RÉQUISITION D'IMMUNOHISTOCHEMIE IMMUNOHISTOCHEMISTRY REQUISITION

Pathologiste demandeur / Requesting Pathologist	
*Nom & Prénom / Last & First Name :	
*No. de Permis / Licence # :	
*Institution :	
*Adresse / Address :	
*No. de téléphone / Telephone No. :	
*Signature :	*Date (AAYY/MM/JD) :

* Renseignements Obligatoires / * Mandatory Information

Information du patient / Patient Information	
* Nom / Last Name :	
* Prénom / First Name :	
* RAMQ :	
Adresse / Address :	
* Date de naissance / Date of birth (AAYY/MM/JD) :	

IMPORTANT

Si possible envoyer une lame de plus au cas où un problème survient.
Envoyer seulement des lames. Les demandes de blocs à couper ne sont pas acceptées.

Double identification des lames envoyées.

Toute demande sans identification conforme sera retournée à vos frais.

Il incombe à l'institution requérante d'ajouter un contrôle sur la lame du patient. Dans le cas contraire, un contrôle de lot ne sera fourni que sur demande.

Pour le tissu de contrôle recommandé sur la lame, se référer au Guide des contrôles positifs pour l'IHC sur le site web du CUSM (<https://muhc.ca/muhc-clinical-laboratories/muhc-laboratory-requisitions-and-consultations>).

If possible please send an extra slide should there be an issue.

Only send slides. Request for blocks to be cut are not accepted.

Double identify the slides sent.

All requests without proper identification will be returned at your expense.

It is the responsibility of the requesting institution to add an on slide control to the patient slide. Otherwise, a batch control will only be

provided if requested. For the recommended on slide control tissue, refer

to the IHC Positive Control Guide in the MUHC website

(<https://muhc.ca/muhc-clinical-laboratories/muhc-laboratory-requisitions-and-consultations>).

No. de pathologie:

Pathology No. :

Nbre de lame(s):

No. of Slide(s):

Spécimen:

Specimen :

BIOMARQUEURS-PARAFFINE / BIOMARKERS -PARAFFIN

☐ MLH1	☐ HER2/NEU	☐ PD-L1 (SP263)
☐ MSH2	☐ ER	
☐ MSH6	☐ PR	
☐ PMS2	☐ BRAF (V600E)	
☐ ROS-1 (SP384)	☐ ALK-Lung (D5F3)	☐ Claudin 18

Veuillez noter que PD-L1 (SP263) est le seul clone disponible au Glen et qu'il n'est réalisé et validé que pour les échantillons de poumon. L'interprétation sera faite par le pathologiste du Glen. Les demandes pour d'autres clones et d'autres indications seront renvoyées au demandeur. /

Note that PD-L1 (SP263) is the only clone available at the Glen and it is only performed and validated for lung specimens. Interpretation will be done by the Glen pathologist. Requests for other clones and other indications will be sent back to the requestor.

Veuillez envoyer 2 lames non colorées et une lame HE. L'interprétation sera faite par le pathologiste de Glen. / Please send 2 unstained slides and an HE slide. The interpretation will be done by the Glen pathologist.

PARAFFINE / PARAFFIN

☐ 2SC	☐ BRG-1 (SMARCA4)
☐ ACTH	☐ C-kit (CD117)
☐ AFP	☐ C-MYC
☐ ALK-HEMA (D5F3)	☐ C4d
☐ Annexin	☐ CAIX
☐ AMACR	☐ Calcitonin
☐ APP	☐ Caldesmon
☐ ARID1A	☐ Calponin
☐ ARID1B	☐ Calretinin
☐ Arginase	☐ Calretinin-Hirschsprung Disease
☐ ATRX	☐ CDX2
☐ B-AMYLOID	☐ CD1a
☐ B-Catenin	☐ CD2
☐ BAP-1	☐ CD3
☐ BCL-2	☐ CD4
☐ BCL-6	☐ CD5
☐ BCL-10	☐ CD7

☐ CD8	☐ CD79A
☐ CD10	☐ CD99
☐ CD15	☐ CD138
☐ CD20	☐ CD163
☐ CD21	☐ Chromogranin
☐ CD23	☐ CKAE1/AE3
☐ CD25	☐ CK5/6
☐ CD30	☐ CK7
☐ CD31	☐ CK8/18
☐ CD34	☐ CK19
☐ CD43	☐ CK20
☐ CD45 (LCA)	☐ Claudin4
☐ CD56	☐ CMV
☐ CD57	☐ Cyclin-D1
☐ CD68 (KP1)	☐ D2-40
☐ CD68 (PGM1)	☐ Desmin
☐ CD71	☐ DOG 1

Information du patient:
Patient Information :

Institution:

☐ E-Cadherin
☐ EMA
☐ ERG
☐ Fumarate Hydratase
☐ FVIII
☐ GATA-3
☐ Gastrin
☐ GFAP
☐ GLUT1
☐ Glutamine Synthetase
☐ Glypican-3
☐ HCG
☐ HCL (DBA.44)
☐ Helico
☐ HEP.B CORE
☐ HEP.B SURFACE
☐ Hepatocyte
☐ HGH
☐ HHV-8
☐ HMB45
☐ HNF1B
☐ HSV I,II Cocktail
☐ IDH1

☐ IgG
☐ IgG4
☐ IgM
☐ Inhibin
☐ INI-1 (BAF47)
☐ Ki-67
☐ Langerin
☐ L-FABP
☐ LH
☐ Lysozyme
☐ Mammaglobin/BRST2 Cocktail
☐ MART1/MELAN A
☐ Merkel Cell
☐ MDM2
☐ MOC-31
☐ MPO
☐ MTAP
☐ MUC4
☐ MUM-1
☐ Myogenin
☐ Napsin A
☐ NEU-N
☐ Neurofilament

☐ NKX3.1
☐ NUT
☐ OCT3/4
☐ OLIG-2
☐ P16ink4
☐ P40
☐ P53
☐ P57KIP2
☐ P63
☐ PAX5
☐ PAX8
☐ Perforin
☐ PIT-1
☐ Polyoma
☐ PRAME
☐ Prolactin
☐ Prostein
☐ PSMA
☐ S-100
☐ SDHB
☐ Serum Amyloid A
☐ SF1
☐ SMA (Actin)

☐ SMAD4
☐ SMI-32
☐ SMM (Myosin)
☐ SOX-10
☐ SOX-11
☐ Spirochete
☐ SSTR2
☐ STAT6
☐ Synaptophysin
☐ TAU (AT8)
☐ TDP-43
☐ TdT
☐ TFE3
☐ Thyroglobulin
☐ Toxoplasma
☐ TPIT
☐ Trypsin
☐ Tryptase
☐ TTF 1
☐ Vimentin
☐ WT-1N

PROBES ISH

ISH PROBES

☐ ISH-EBER Probe
☐ ISH-Kappa Probe
☐ ISH-Lambda Probe

CHROMOGÈNE ROUGE

RED CHROMOGEN

☐ A-SYNUCLEIN (RED)
☐ MART1/MELAN A (RED)
☐ S-100 (RED)
☐ HMB45 (RED)
☐ SOX-10 (RED)

DOUBLE COLORATIONS

DUAL STAINS (BROWN + RED)

☐ CD31 + CKAE1/AE3
☐ D2-40 + CKAE1/AE3
☐ Ki-67 + Mart-Mel
☐ P40 + CK5/6
☐ PIN3-DUAL
☐ TTF1 + Napsin A

CONGÈLATIONS / FROZEN SECTIONS

BILAN DYSTROPHIN PANEL

DYSTROPHIN PANEL

☐ Caveolin 3
☐ A-Dystroglycan
☐ B-Dystroglycan
☐ Dystrophin-1
☐ Dystrophin-2
☐ Dystrophin-3
☐ Dysferlin
☐ Merosin
☐ A-Sarcoglycan
☐ B-Sarcoglycan
☐ D-Sarcoglycan
☐ G-Sarcoglycan
☐ Spectrin
☐ Utrophin

BILAN MUSCULAIRE

MUSCLE PANEL

☐ MAC
☐ Myosin HC Fast / Laminin (DUAL STAIN)
☐ Myosin HC Slow / Laminin (DUAL STAIN)
☐ Myosin HC Neonatal / Laminin (DUAL STAIN)

BILAN IMMUNOFLUORESCENCE

IMMUNOFLUORESCENCE PANEL

☐ FITC IgA
☐ FITC IgG
☐ FITC IgM
☐ FITC C1q
☐ FITC C3
☐ FITC Kappa
☐ FITC Lambda
☐ IF C4D (Indirect)